

# Session 3: Beginning manual curation Day 2


Genome Reference Informatics Team (GRIT) Wellcome Sanger Institute - Tree of Life





- Analysis pipelines
- Treeval
- Curation pretext pipeline
- Curation tools
  - Jbrowse

  - Rapid curation workflow (tomorrow)
    How to produce a curated fasta file
- How to use PretextView





# TreeVal as an aid for curation


**Jbrowse** 

# TreeVal Browser




• All projects coming into the curation team now go through the TreeVal Pipeline (github.com/sanger-tol/treeval)

#### All draft assemblies



Might not be enough for curation

# **Evidence data for curation** e.g HiC maps, Coverage, telomere, repeat tracks, gene alignments ranka ara dalaman ara mana manika dalaman manika






TreeVal Browser

In depth information is provided

## TreeVal Browser





# Why Jbrowse?


- Jbrowse provides the perfect platform for visualisation of assemblies:
- Simple to use, web-based, lightweight, interactive and scalable
- It's an extensible platform that utilizes multiple data types in a variety of genomic data formats such as FASTA, VCF, BAM, BED/BIGBED, HIC, CRAM, PAF and GFF3
- Highly customisable

How to start using TreeVal (Jbrowse)?

http://jbrowse.tol.sanger.ac.uk



#### Type your TOLID here





### How to start using TreeVal?




| TreeVal    |           |             |
|------------|-----------|-------------|
| ToLID \$   | jBrowse   |             |
| bmerAlb1_1 | ~         |             |
| bMerAlb1_1 | Available | <u>View</u> |

Showing 1 to 1 of 1 rows


#### How to start using TreeVal?

alignment




# Which data types do we align to our assemblies?





# Gene alignment data






# Synteny data









# How to use PretextView





Download PretextView from: https://github.com/sanger-tol/PretextView





# How to use PretextView



- Click 'Load Map' and open a map
- Keyboard Shortcuts
- Press 'U' to open menu
- Press 'E' to enter the Edit mode and make edits
  - Press 'Q' to undo an edit
  - Press 'W' to redo an edit
  - Click middle button to flip (invert) a scaffold

#### Press 'S' to enter the Scaffold mode: Paint scaffolds blocks as chroms

Click middle button to delete painting

Press 'A' to paint a block of scaffolds at once

Press Shift + 'D' to delete all the painting in the map

# How to use PretextView



- Press 'M' to enter Metadata Tag mode: add tags such as haplotigs, unloc, HAP1, HAP2, X, Z, Contaminant...
  - Use arrow keys in your keyboard to select the active tag
  - Press Shift + 'D' to delete all tags in your map
- Press 'I' to show ID bars: delimit and show orientation of the scaffolds
- Press 'G' to switch the grids on and off
- Press 'B' to show scale bars at the top
- Press 'T' to switch the tool tip on and off: show coordinates in the map and info about active tracks
- Arrow keys left/right: adjust contrast
- Arrow keys up/down: change color scheme in the map
- Press 'W' to enter Waypoint edit mode: mark specific regions